



- Slice Anchor Kits are available for all Series 20 tissue slice recording chambers.
- Kits for RC 26, RC-26G, RC 27, and RC-27L, (large bath chambers) contain three stainless steel anchors with 0.1 mm thick lycra threads at 1 mm, 1.5 mm, and 2.0 mm spacing.
- Kits for the RC -22 and RC- 22C (smaller slice chambers) contain two stainless steel anchors with lycra threads at 1 mm and 1.5 mm spacing.
- Stainless steel material allows for anchors to be adjusted for a tighter or easier fit into chamber
- Each kit comes with nylon mesh in two sizes (112 micron and 160 micron) that may be cut and used in conjunction with the slice anchor.

**Ordering information**

SHD-22KIT	Slice anchor kit for RC-22 chamber
SHD-22CKIT	Slice anchor kit for RC-22C chamber
SHD-26KIT	Slice anchor kit for RC-26 chamber
SHD-26GKIT	Slice anchor kit for RC-26G chamber
SHD-27KIT	Slice anchor kit for RC-27 chamber
SHD-27LKIT	Slice anchor kit for RC-27L and RC-29 chambers

